

Cournon d'Auvergne, June 05<sup>th</sup>, 2020

# easyret®

## WHEN FIBER LASER TECHNOLOGY MEETS IMAGING




Dear Partners,

We are pleased to inform you that **Mosar®**, our new **High Definition Imaging System** for Easyret is now available for order.


Composed of a HD camera, a second screen and a powerful computer, Mosar® is an option for Easyret® featuring:

1. A **co-observation teaching mode** for live viewing of laser procedures.
2. An **advanced mode** allowing the operator to:
  - Import diagnosis images facilitating the laser treatment planning
  - Prepare, print and record advanced treatment reports including fundus and diagnosis images
  - Take pictures or record treatment videos for presentation and training purposes
3. A **library mode**, to manage all the generated images, videos and treatment reports which can be exported on a USB drive or a local network.


### 1. CO-OBSERVATION MODE



### 2. ADVANCED MODE




### 3. LIBRARY MODE



FOLLOW US!



# EASYRET SOFTWARE 2.01



Along with the introduction of the Mosar® system, we are also happy to announce the availability of a new software release for Easyret® (version 2.01).

This new software featuring our new advanced titration procedure designed by Victor Chong, replaces the previous 2.0 version and is now compatible with all the Easyret laser versions:

- It is provided as a standard with all new Easyret® laser orders.
- All the Easyret® lasers (with our without HD screen) already installed are fully compatible and eligible for software update.

Feel free to contact our service department such you have any question.

## AVAILABILITY

Due to the coronavirus pandemic and its impact on our manufacturing capacities Mosar® will be available in limited quantities. The production of our new imaging system should return to its normal level by **September 2020**.

As result, we inform you that during the first 3 months of its commercialization, Mosar® will only be available for order in combination with a new Easyret® system:

### Easyret laser order codes:

- EASYRET® Laser HAAG STREIT TYPE: **PABX0041HSSI**
- EASYRET® Laser HAAG STREIT TYPE with LIO port: **PABX0041HSDI**
- EASYRET® Laser ZEISS TYPE: **PABX0041ZESI**
- EASYRET® Laser with LIO port ZEISS TYPE: **PABX0041ZEDI**
- **Mosar® order code: PAEX0041MOS**

But don't worry, as announced originally, Mosar will be compatible and upgradable with all the already installed Easyret lasers. Additional information will be announced in **September 2020**.

### FOLLOW US!





# MARKETING TOOLS

For your information, a marketing kit is now available on our extranet webpage including:

1. A sales brochure in English, Spanish and Polish (Mosar® system, new clinical content etc.)
2. A promotional in English, Spanish and Polish
3. A roll up in English
4. An English advertisement

1.



2.



3.



4.



We would appreciate if you could make sure to share this information with your customers. If you need any further details, please contact your area sales manager, or your sales assistant.

We wish you a lot of success and happy sales!

The Marketing Team

**Internal communication and confidentiality Notice:** This document is confidential and contains proprietary information and intellectual property of Quantel Medical. Please be aware that copying or distribution of this document and the information contained therein is **strictly prohibited outside your company.**

